

Hei-SHAKE Reciprocating Core

Technical datasheet



| | |
|---|---|
| Dimensions (W × H × D) | 356 x 129 x 403 mm |
| Weight | approx. 14 kg |
| Maximum permissible load | 15 kg |
| Drive | EC Motor |
| Speed range | 1 - 300 rpm |
| Acoustic pressure | < 85 dB (A), in accordance with IEC 61010 |
| Rated voltage | 1/N/PE 100-240 V, AC 50/60 Hz |
| Overvoltage category | N/A |
| Degree of pollution | N/A |
| Power input | 90 W |
| Operating temperature | 5 °C – 31 °C at up to 80 % rel. humidity |
| Installation altitude | up to 2,000 m asl |
| Overheat protection | N/A |
| Speed control | N/A |
| Usable area | (W x D) 294 × 263 mm |
| Movement | reciprocating (horizontal shaking motion) |
| Orbit | N/A |
| Fuse | N/A |
| Permissible mains voltage fluctuations | N/A |
| Bottle dimensions | N/A |
| Interfaces | N/A |

Protection Class

IP32 (acc. to EN 60529)