

Reax 20 - 8 (1~16rpm)
Technical datasheet



Dimensions (W × H × D)	770 × 520 × 465
Weight	28 kg
Maximum permissible load	30 kg
Drive	Asynchronous motor
Speed range	1 – 16 rpm; Speed display: scale
Protection class (EN 60529)	IP21
Acoustic pressure	< 45 dB (A), in accordance with IEC 61010
Rated voltage	230 V, 50/60 Hz 115 V, 60 Hz
Protection class	I
Overvoltage category	n.a.
Degree of pollution	n.a.
Power input	280 W
Operating temperature	5 °C – 31 °C at up to 80 % rel. humidity 32 °C – 40 °C at up to 50 % rel. humidity (decreasing linearly)
Installation altitude	up to 2,000 m asl
Overheat protection	self-resetting
Speed control	electronic
Usable area	n.a.
Movement	vertically rotating
Orbit	n.a.

Fuse

T3, 15 A (230 V) | T6, 25 A (115 V)

**Permissible mains voltage
fluctuations**

n.a.

Bottle dimensionsMaximum bottle diameter: 136 mm | Maximum bottle
height: 270 mm**Interfaces**

n.a.