

Geno™ II

Technical datasheet



Dispense Flow Rate	3 - 7 L/min
Water Quality (at 25° C)	1 - 15 MΩ·cm
Pressure	0.1 - 6 Bar/1.45 - 86 PSI
Conductivity	<1400µS/cm
Temperature	4 - 30°C/39.2 - 86°F
Free Chlorine	<2 ppm
pH Value	4 - 10
Dimensions (H x W x D)	500 x 490 x 290 mm/19.7" x 19.3" x 11.5"
Installed Weight	23 kg/50.7 lbs
Power Supply	100 - 240 v ~; 50/60 Hz; 3.0A (max)
Silt Density Index	< 3 SDI