

## Hei-CHILL Peltier Technical datasheet



<b>Dimensions (W × H × D)</b>	27 x 43 x 43 cm
<b>Weight</b>	15 kg
<b>Heating power</b>	>1000 W
<b>Temperature measurement accuracy</b>	0.05 °C
<b>Acoustic pressure</b>	63 db
<b>Control</b>	Panel: 5" Touchscreen
<b>Working temperature range</b>	-10 - 70 °C
<b>Cooling capacity</b>	600 W
<b>Interfaces</b>	Ethernet, RS 232, RS 485
<b>Flow rate</b>	max. 3 l / min