

## Hei-SHAKE Rock Core 2000 (5°) Technical datasheet



<b>Dimensions (W × H × D)</b>	405 × 219 × 430 mm
<b>Weight</b>	approx. 14 kg
<b>Maximum permissible load</b>	max. 15 kg
<b>Drive</b>	EC motor
<b>Speed range</b>	0 - 80 rpm
<b>Acoustic pressure</b>	< 85 dB (A), in accordance with IEC 61010
<b>Rated voltage</b>	1/N/PE 100-240V AC 50/60 Hz
<b>Protection Class</b>	IP32 (acc. to EN 60529)
<b>Power input</b>	90 W
<b>Usable area</b>	aprox. 390 × 340 mm
<b>Movement</b>	Rocking
<b>Interfaces</b>	1 × RS232 (Sub-D, 9 pole)
<b>Tilt Angle</b>	5°
<b>Operating and Storage Temperature</b>	5 °C – 31 °C at up to 80 % rel. humidity 32 °C – 40 °C at up to 50 % rel. humidity (decreasing linearly)
<b>Installed Altitude</b>	up to 2,000 m asl