

Findenser™ Mini Cône NS24, Douille NS24

Technical datasheet



| | |
|--------------------------------|--------------------|
| Poids | 670 g |
| Length × Diameter | 275 mm × 72 mm |
| Joints | B14, B19, B24 |
| Inner or wetted parts | Borosilicate glass |
| Finned jacket | Anodised aluminium |
| Plug & Seal | Acetal & FKM |
| Thermofluid | Water |
| Max solvent bpt | 155 °C |
| Oven drying | Max Temp. 60 °C |
| Dishwasher | Max Temp. 50 °C |
| Autoclaving | Do Not Autoclave |
| Ambient operating temp. | 0 °C to 24 °C |