

Hei-CHILL 700 Pro

Technical datasheet



Dimensions (l × h × p)	External: 355 x 455 x 505 mm Installation (W x D): 355 x 755 mm
Poids	39 kg
Classe de protection (EN 60529)	IP 20
Niveau de pression acoustique	49 dB(A)
Type de raccordement	Electrical: 230 V / 50 Hz
Alimentation	500 W
Précision de réglage	at -10°C: +/- 0,5 K
Working temperature range	-10°C to +40°C
Display resolution	0,1 K
Cooling capacity	at 20°C: 700 W at 10°C: 550 W at 0°C: 350 W at -10°C: 120 W
Circulation pump	Flow rate [max.] 10,5 l/min Flow pressure [max.] 2,8 bar
Hose connections	3/8" female thread
Coolant tank	5,5 to 9,0 l
Permitted ambient temperature range	5°C to 40°C
Max. temperature of coolant	60°C
Current [max.]	3,5 A
Class. acc. DIN 12876	I / NFL
Max. operating pressure	25 bar
Refrigerant	R290 (GWP 3*)