

Hei-VAC Vario Control Technical datasheet



Dimensions (l × h × p)	236 x 167 x 196 mm
Poids	5.4 kg
Classe de protection (EN 60529)	IP 20
Type de raccordement	Coolant: hose nozzle DN 6-8 mm
Température de service	+10 - +40°C
Protection du moteur	temperature sensor on the pcb (current limitation)
Storage temperature range	-10 - +60°C
Puissance	Electrical: 0.16 kW
Maximum pumping speed	1.7 m ³ /h
Ultimate vacuum (absolute)	without gas ballast: 2 mbar with gas ballast: 4 mbar
Maximum permissible pressure	1.1 bar
Permissible relative atmospheric moisture during operation (no condensation)	30% - 85%
Maximum no-load speed	2400 min ⁻¹
Inlet	hose nozzle DN 8 mm
Outlet	hose nozzle DN 8 mm
Maximum permissible range of supply voltage (±10%)	100-120 V~ 50/60 Hz 200-240 V~ 50/60 Hz