

## Findenser™ Mini NS24 Cone, NS24 Socket

### Technical datasheet



<b>Weight</b>	670 g
<b>Length × Diameter</b>	275 mm × 72 mm
<b>Joints</b>	B14, B19, B24
<b>Inner or wetted parts</b>	Borosilicate glass
<b>Finned jacket</b>	Anodised aluminium
<b>Plug &amp; Seal</b>	Acetal & FKM
<b>Thermofluid</b>	Water
<b>Max solvent bpt</b>	155 °C
<b>Oven drying</b>	Max Temp. 60 °C
<b>Dishwasher</b>	Max Temp. 50 °C
<b>Autoclaving</b>	Do Not Autoclave
<b>Ambient operating temp.</b>	0 °C to 24 °C