

# Hei-VOLUME Distimatic Pro Industrial with adaptor for A-series glassware

## Technical datasheet



|                                    |  |
|------------------------------------|--|
| <b>Dimensions (W × H × D)</b>      | 186 × 429 × 521 mm, including control panel  |
| <b>Weight</b>                      | approx. 17 kg  |
| <b>Protection class (EN 60529)</b> | IP42   |
| <b>Rated voltage</b>               | 100 – 240 V AC; 50/60 Hz   |
| <b>Connection</b>                  | L+N+PE   |
| <b>Protection class</b>            | I  |
| <b>Overvoltage category</b>        | II   |
| <b>Degree of pollution</b>         | 2  |
| <b>Power input</b>                 | 1,500 W max  |
| <b>Operating temperature</b>       | 5 °C – 31 °C at up to 80 % rel. humidity   32 °C – 40 °C at up to 50 % rel. humidity (decreasing linearly) |
| <b>Installation altitude</b>       | up to 2,000 m asl  |
| <b>Fuse</b>                        | 15 A   |