

Hei-SHAKE Salto Core Technical datasheet



| | |
|---|--|
| Dimensions (W × H × D) | 1.110 × 520 × 465 mm |
| Weight | approx. 31 kg |
| Maximum permissible load | 30 kg |
| Drive | Gear stepper motor |
| Speed range | 1 – 32 rpm (clockwise/counter-clockwise rotation) |
| Acoustic pressure | < 60 dB (A), in accordance with IEC 61010 |
| Rated voltage | 1/N/PE 100-240V AC 50/60 Hz |
| Overvoltage category | II |
| Degree of pollution | 2 |
| Power input | max. 90 W |
| Operating temperature | 5 °C – 31 °C at up to 80 % rel. humidity 32 °C – 40 °C at up to 50 % rel. humidity (decreasing linearly) |
| Installation altitude | up to 2,000 m asl |
| Overheat protection | n.a. |
| Speed control | n.a. |
| Usable area | n.a. |
| Movement | Overhead shaker |
| Orbit | n.a. |
| Fuse | n.a. |
| Permissible mains voltage fluctuations | n.a. |

Bottle dimensions

n.a.

Interfaces

1 × RS232 (Sub-D, 9 pole)

Protection Class