

Hei-VAC Vario Control

Technical datasheet



| | |
|---|---|
| Dimensions (W × H × D) | 236 x 167 x 196 mm |
| Weight | 5.4 kg |
| Protection class (EN 60529) | IP 20 |
| Connection | Coolant: hose nozzle DN 6-8 mm |
| Operating temperature | +10 - +40°C |
| Motor protection | temperature sensor on the pcb (current limitation) |
| Storage temperature range | -10 - +60°C |
| Power | Electrical: 0.16 kW |
| Maximum pumping speed | 1.7 m ³ /h |
| Ultimate vacuum (absolute) | without gas ballast: 2 mbar with gas ballast: 4 mbar |
| Maximum permissible pressure | 1.1 bar |
| Permissible relative atmospheric moisture during operation (no condensation) | 30% - 85% |
| Maximum no-load speed | 2400 min ⁻¹ |
| Inlet | hose nozzle DN 8 mm |
| Outlet | hose nozzle DN 8 mm |
| Maximum permissible range of supply voltage (±10%) | 100-120 V~ 50/60 Hz 200-240 V~ 50/60 Hz |