

## Hei-VAC Vario Station Technical datasheet



<b>Dimensions (W × H × D)</b>	294 x 192 x 299 mm
<b>Weight</b>	6 kg
<b>Protection class (EN 60529)</b>	IP 20
<b>Connection</b>	n.a.
<b>Operating temperature</b>	+10 - +40°C
<b>Motor protection</b>	temperature sensor on the pcb (current limitation)
<b>Storage temperature range</b>	-10 - +60°C
<b>Power</b>	Rated: 0.16 kW
<b>Maximum pumping speed</b>	1.7 m <sup>3</sup> /h
<b>Ultimate vacuum (absolute)</b>	without gas ballast: 5 mbar   with gas ballast: 7 mbar
<b>Maximum permissible pressure</b>	1.1 bar
<b>Permissible relative atmospheric moisture during operation (no condensation)</b>	30% - 85%
<b>Maximum no-load speed</b>	0 - 2400 min <sup>-1</sup>
<b>Inlet</b>	hose nozzle DN 8/10 mm
<b>Outlet</b>	hose nozzle DN 8 mm
<b>Maximum permissible range of supply voltage (±10%)</b>	100-120 V~ 50/60 Hz   200-240 V~ 50/60 Hz