

Hei-VOLUME Distimatic Pro 24/7 Industrial con adaptador para vidrio R

Technical datasheet



| | |
|------------------------------------|--|
| Dimensions (W × H × D) | 186 × 429 × 521 mm, including control panel |
| Weight | approx. 17 kg |
| Protection class (EN 60529) | IP42 |
| Rated voltage | 100 – 240 V AC; 50/60 Hz |
| Connection | L+N+PE |
| Protection class | I |
| Overvoltage category | II |
| Degree of pollution | 2 |
| Power input | 1,500 W max |
| Operating temperature | 5 °C – 31 °C at up to 80 % rel. humidity 32 °C – 40 °C at up to 50 % rel. humidity (decreasing linearly) |
| Installation altitude | up to 2,000 m asl |
| Fuse | 15 A |